

In-situ Grazing Incidence Small-Angle-X-ray-Scattering Studies Using High Brilliance Microfocus Sources

Jörg Wiesmann¹, Andreas Stricker¹, Andreas Kleine¹, Martin Hodas², Karol Vegso², Peter Siffalovic²

¹Incoatec GmbH, Max-Planck-Str. 2, 21502 Geesthacht, Germany; Email: sales@incoatec.de; ²Institute of Physics, Slovak Academy of Sciences, Dubravská cesta 9, 845 11 Bratislava, Slovakia

Grazing Incidence Small Angle Scattering (GISAXS)

Measures diffuse scattering in reflection geometry at grazing incidence (α_i) and grazing reflectance (α_r) angles

Diffuse scattering contains information about


- surface and subsurface structure
- island dimensions and 3-D arrangements
- roughness
- pore diameter

Advantages:

- non-destructive, allowing inspection of buried clusters and interfaces
- in contrast to microscopy average over sample surface/volume
- working under all conditions (enables in-situ screening/control)
- changing the incidence angle controls the penetration depth

Disadvantages:

- transformation from reciprocal to real space needs an appropriate model
- we always measure statistical averages over ensemble of entities





The potential of our microfocus source μS in GISAXS studies in the home-lab is demonstrated in an overview of representative experimental setups and results. The studies take advantage of the brilliance and outstanding beam quality of the low-maintenance μS . It is shown how the μS can be used to achieve excellent results in the investigation of in-situ thin film deposition in UHV chambers by using GISAXS or of the structure of oriented two-dimensional liquid crystalline samples.


In-situ GISAXS during thin film growth

By using in-situ GISAXS in the home-lab we investigated how multilayers grew during thin film deposition. This kind of experiments is typically done only at synchrotrons. With an μS it is now also feasible in the home-lab.


μS for in-situ GISAXS with Pilatus 200K




The Dual Ion-beam sputtering unit in Bratislava was upgraded with an in-situ GISAXS set-up. As a source an μS with a special collimating optics for SAXS is mounted on a Hexapod. Together with the 2-dim detector Dectris Pilatus 200K dynamic measurements during thin film growth become feasible.




Reciprocal space map of 10 periods W/B₄C multilayer mirror with 1.5 nm period thickness measured ex-situ by GISAXS in deposition chamber



Time resolved evolution GISAXS reciprocal space map of the 10x W/B₄C multilayer mirror with visible Bragg peak and Kiessig fringes



GISAXS reciprocal space map of the 40 x Mo/Si multilayer mirror with period 6.9 nm




Time resolved evolution GISAXS reciprocal space map of the 40x Mo/Si multilayer mirror with visible Bragg peaks

GISAXS plots show the perfect growth of the multilayers. Even thin films with a total thickness in the range of 15 nm could be measured. The time resolved evolution of the specular signal enables the measurement of the Bragg peaks and the Kiessig Fringes dynamically.

NEW and in progress: In-situ time-resolved GISAXS of metal films on graphene

- This method revealed kinetics of Cu cluster growth on CVD graphene.
- It allows rapid optimization of metal deposition processes in laboratory conditions.
- Further growth studies of Au, Ag,..., on graphene surface are in progress.



Upgrading Existing Diffractometers with the Incoatec Microfocus Source μS

You have a Bruker AXS, Marresearch, Nonius, Rigaku, Huber or some other system? Incoatec offers a unique possibility to upgrade your existing diffractometer by installing the high-performance, air-cooled and low-power microfocus source μS .

Your upgrade benefits:

- No maintenance, only single phase power and no water cooling required
- 3 years warranty
- Maximum installation down time of only 2-4 days
- New safety concept development on request
- Full compliance with European Machinery Directive 2006/42/EC

Your upgrade options:

- Source, optics and scatterless slits
- Single source upgrade for XRD, XRR, (GI)SAXS and many more applications
- Dual wavelength setup by adding μS as complementary source
- Cu, Mo, Ag, Co and Cr radiation (others on request)



The three generations: μS , $\mu\text{S}^{\text{high Brilliance}}$ and μS 3.0



μS + SCATEX upgrade on a customized SAXS setup in Hamburg



Replacement of Rigaku RU-200 generator in Boulder, USA




Huber system for SAXS in Tamkang, Taiwan

In-situ GISAXS of Nano-particles on liquid surfaces


GISAXS measurements were performed with the μS /Dectris system at the Institute of Physics at the Slovak Academy of Science in Bratislava. Silver particles on a Langmuir film were analyzed at different surface pressures which were applied by means of a reduction of the surface area.

μS with Dectris pixel detector

For rapid grazing incidence small-angle-scattering measurements (GISAXS) of liquid samples our μS equipped with an optics with 5mrad divergence was combined with a Dectris Pilatus detector. The set-up including the alignment can be changed from liquid to capillary samples in less than 30 min.




Analyzing the particle diameter: SAXS of diluted solution of Ag nanoparticles. Result: Metal core, radius = 2.9 ± 0.3 nm




Experimental for GISAXS:

- angle of incidence: 0.2 deg
- 180 s measurement time
- aperture 350 μm
- the surface was pressed with 0 up to 26 mN/m



SAXS patterns at three different surface pressures:

- 0 mN/m: unpressed
- 16 mN/m: intensity increases
- 26 mN/m: crystal-like peaks appear



- Unpressed surface: islands of nanoparticles are swimming on the surface without connection
- Increasing surface pressure: islands coalescence
- vertical formation of hexagonal layers

Ordering phenomena could be observed in-situ during a change of surface pressure. The particles were transformed from single islands to an almost vertically ordered structure of connecting particles.